

912-0007

YC-100™ Test Kit (diploid mode)Certificate & Package Insert for lot **0616-20**. Expire by **Sep 2017**

	ChemoMetec Master Instrument	Accept Criteria	
		Average of 5 analysis	CV%
BLANK (ZERO)	0.00 x 10 ⁶ beads/ml	0.00 – 0.20 x 10 ⁶ beads/ml	
Fluorescent beads	0.71 x 10 ⁶ beads/ml	0.51 – 0.91 x 10 ⁶ beads/ml	Max. 8%

This lot has been tested according to the internal ChemoMetec quality procedures. The lot meets or exceeds all product specifications.

Certificate released by Cecilia Klyver/QC Manger, in June 2016.

Contents	One vial containing 1.0 ml fluorescent bead suspension and one vial containing 1.5 ml of a bead free suspension (ZERO). The vials contain sufficient volume for at least five analyses. Discard the remaining suspension after use or when expiry date is reached.
Application	The YC-100™ Test Kit is intended to be used for routine check of the NucleoCounter® YC-100™ instrument. It is recommended to perform a test of the NucleoCounter® twice per year. The beads of the kit are fluorescent and will appear like yeast cells in the image presented by the NucleoView™ software.
Procedure	<ul style="list-style-type: none"> - The NucleoCounter® YC-100™ must be in "diploid" mode when using the YC-100™ Test Kit. - Shake the Fluorescent bead vial vigorously for 5 seconds before use. In order to obtain a single-bead suspension it is extremely important to shake the vial sufficiently prior to use. Vortex mixing CANNOT replace shaking! - After shaking, let the vial stand for 3-5 minutes to allow the foam to settle. - BLANK: Remove the lid from vial, load a NucleoCassette™ with the solution and then carry out an analysis on the NucleoCounter®. Repeat the test procedure five (5) times using five NucleoCassettes™. - Fluorescent beads: Invert the vial gently, remove the lid, load a NucleoCassette™ with the Fluorescent bead suspension and then carry out an analysis on the NucleoCounter®. Repeat the test procedure five (5) times using five NucleoCassettes™. Close the vial and invert it gently before each measurement in order to avoid settling of the beads. - Calculate the average count in millions/ml for each suspension. Discard up to one outlier for each calculation of average counts. - If the average counts are within the specified range (refer to the top of this document) the NucleoCounter® YC-100™ is working satisfactory. If the average counts are not within the acceptable range, this can be due to problems with the instrument, the NucleoCassettes™, the Test Kit or the handling of the system. Pay special attention to shaking of the bottles prior to use.
Stability	The storage stability after production is at least 9 months if stored at 2-6°C. See expiry dates on the kit labels and on top of certificate.
Storage	Store the kit at 2-6°C until opening. DO NOT FREEZE! After opening the kit shall be used immediately. After use discard the kit.
Shipment	The kit can be shipped at ambient temperature (0-50°C for max. 10 days).
Safety Information	This material is not considered hazardous. However, exercise due care and wash hands after usage as a precaution. Do not eat or drink the product.